a *Johnson & Johnson* company

# **ABSOLUTE™ ABSORBABLE INTERFERENCE SCREW SYSTEM**

**DePuy Mitek, Inc.**  
325 Paramount Drive  
Raynham, MA 02767, USA  
1-800-356-4835

**CE 0086**

EC	REP
----	-----

**DePuy International Ltd.**  
St Anthony's Road  
Leeds LS11 8DT England  
Tel: +44 113 270 0461  
Fax: +44 113 272 4101

**P/N: 105494      Rev: H      02/08**

# ENGLISH

## ABSOLUTE Absorbable Interference Screw System

### Instructions for Use

#### DESCRIPTION

The DePuy Mitek ABSOLUTE Absorbable Interference Screw is an absorbable, tapered, threaded fastener for use in interference fixation of soft tissue grafts or bone-tendon grafts.

The ABSOLUTE Absorbable Interference Screw is packaged one per pouch.

#### INDICATIONS

The DePuy Mitek ABSOLUTE Absorbable Interference Screw is intended for attachment of soft tissue grafts or bone-tendon grafts to the tibia and/or femur during cruciate ligament reconstruction surgeries.

#### CONTRAINdications

1. This device is not cleared for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
2. Insufficient quantity or quality of bone, comminuted bone surfaces, or pathologic bone conditions such as cystic change or severe osteopenia that would impair the ability of the ABSOLUTE Absorbable Interference Screw to securely fixate to the bone.
3. Physical conditions which would eliminate, or tend to eliminate, adequate implant support or retard healing, i.e., blood supply limitations, infections, etc.
4. Implant/material sensitivity/foreign body sensitivity.
5. Conditions which tend to pre-empt the patient's ability or healing period, such as senility, mental illness or alcoholism.

## **WARNINGS**

1. Inspect instruments for damage prior to use. Replace damaged or worn instruments. Do not attempt to repair.
2. If resistance is felt during the insertion of an ABSOLUTE Absorbable Interference Screw over a guide wire, stop and confirm that the guide wire is not entrapped. If entrapped, back out the screw and withdraw the guide wire.

## **PRECAUTIONS**

1. The DePuy Mitek ABSOLUTE Absorbable Interference Screw should not be used without prior review of the procedure and Instructions for Use.
2. The DePuy Mitek ABSOLUTE Absorbable Interference Screw implant is supplied STERILE and is intended for single use only. Do not resterilize.
3. It is mandatory to completely seat the ABSOLUTE Interference Screw onto the correct driver for proper implant delivery.

## **CONTENTS**

Contents are STERILE unless packaging is damaged or opened.

## **ADVERSE EFFECTS**

Adverse effects of absorbable implanted devices include mild inflammatory and foreign body reactions.

## **INSTRUCTIONS FOR USE**

### **(*for femoral fixation*)**

1. With the knee hyperflexed beyond 90°, insert a 1.1mm Nitinol guide wire through the anterior medial portal or patellar defect, into the femoral tunnel beside the graft.
2. Place the appropriate size Interference Screw on the tip of the Cannulated Hex Driver. Advance the Hex Driver over the Nitinol guide wire into the gap between the graft and the tunnel wall.
3. Rotate the Hex Driver clockwise until it is counter sunk slightly inside the tunnel. Remove the Hex Driver and guide wire. If the nitinol guide wire becomes entrapped, rotate the Hex Driver counter clockwise to back out the screw until the guide wire straightens. Pull back on the guide wire, and readvance the screw.

**(for tibial fixation)**

1. With tension on the graft, insert a guide wire into the gap between the graft and the tibial tunnel wall.
2. Follow steps 2 & 3 above.

**CLEANING AND STERILIZATION INSTRUCTIONS**

Cleaning and sterilization instructions (reference IFU #108410) are provided with the instrument and must be reviewed prior to use.

**STORAGE**

The boxes should be stored under cool and dry conditions.

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Patents pending. All rights reserved.

- A** JOHNSON & JOHNSON Medical GmbH  
DePuy  
Gunoldstrasse 16  
1190 Wien
- B** JOHNSON & JOHNSON Medical BV  
Eikelenbergstraat 20  
B-1700 Dilbeek
- CH** JOHNSON & JOHNSON Medical  
Rotzenbuehlstrasse 55  
CH-8957 Spreitenbach
- CZ** JOHNSON & JOHNSON s. r. o.  
Karla Englise 3201/6  
150 00 Parha 5  
Česká republika
- D** ETHICON GmbH  
Mitek Division  
Oststrasse 1  
D-22844 Norderstedt
- DK** JOHNSON & JOHNSON  
Bregnerødvej 133 - 3460 Birkerod
- E** JOHNSON & JOHNSON  
Productos Profesionales  
Div. ETHICON/Mitek  
Paseo de las Doce  
Estrellas, 5-7  
Campo de las Naciones  
28042 Madrid
- F** ETHICON S.A.S.  
TSA 81002  
1, rue Camille Desmoulins  
92787 Issy les Moulineaux CEDEX 9
- FIN** JOHNSON & JOHNSON  
Metsanneidonkuja 10 C2130 ESPOO
- GB** ETHICON Ltd.  
Bankhead Ave., Edinburgh, EH11 4HE

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**RU**

JOHNSON & JOHNSON LLC  
43, Bolshaya Tulskaya Street  
Moscow, 115191  
Russia

**SK**

JOHNSON & JOHNSON Slovakia, s.r.o.  
Jakubovo námestie 13  
811 09 Bratislava  
Slovenská republika

**S**

JOHNSON & JOHNSON AB  
Staffans Vag 2  
S-191 84 Sollentuna

**TR**

Ethicon Limited  
Turkiye Irtibat Office (J&J Medical)  
Cumhuriyet Cd.  
Gursel Plaza No. 99  
Kat: 5 81610 Kavacik

**CN**

JOHNSON & JOHNSON Medical (Shanghai) Ltd.  
4th. Floor, Cheng Kai International Tower,  
355 Hong Qiao Road,  
Shanghai 200030, China

**TW**

JOHNSON & JOHNSON (HK) Ltd.  
Taiwan Branch  
6F, 319, Sec. 2, Tun Hwa South Road  
Taipei 106

**KR**

JOHNSON & JOHNSON Medical Korea Ltd.  
Medical Korea Ltd.  
14F kukje Building  
191 Hangangro-2Ga  
Yangsan-Gu  
Seoul, Korea

Australia: JOHNSON & JOHNSON Medical Pty. Ltd.  
1-5 Khartoum Road  
North Ryde, NSW 2113

Canada: JOHNSON & JOHNSON Medical Products  
200 Whitehall Drive  
Markham, Ontario L3R 0T5

S. Africa: JOHNSON & JOHNSON (Pty) Ltd.  
2 Medical Street, (PO Box 273)  
Halfway House  
Midrand Gauteng  
S. Africa 1685

**QTY**

Number of Units  
Número de unidades  
Numero di unità  
Anzahl  
Nombre d'unités  
Aantal eenheden  
Número de unidades  
Antal enheder  
Antall enheter  
Yksikköjen lukumäärä  
Αριθμός μονάδων  
Antal enheter  
Počet kusů  
Počet kusov  
Ilość  
Число единиц  
Egységek száma  
Birim Adedi  
**件数**  
裝置數量  
수량

Recycle Packaging  
Reciclar el envase  
Riciclare l'imballaggio  
Verpackung zum Recycling verwenden  
Recyclez l'emballage  
Verpakking opnieuw gebruiken  
Reciclar embalagem  
Genbrugsemballage  
Gjenvinnbar emballasje  
Kierrätä pakkaus  
Ανακυκλώστε τη συσκευασία  
Återvinn förpackningen  
Recyklujte tento obal  
Recyklujte tento obal'  
Opakowanie nadające się do recyklingu  
Упаковка подлежит вторичной переработке  
Újrafelhasználandó csomagolás  
Ambalajı geri kazanın  
**对包装进行再循环**  
**對包裝進行再迴圈**  
**포장 재활용 가능**

MADE IN

LOT

STERILE EO

CE 0086