a *Johnson & Johnson* company

# VERSALOK™ Anchor

**DePuy Mitek, Inc.**

325 Paramount Drive

Raynham, MA 02767, USA

1-800-356-4835

**C** **€**<sub>0086</sub>

**DePuy International Ltd.**

St Anthony's Road, Beeston

Leeds LS11 8DT

United Kingdom




Tel: +44 113 270 0461

Fax: +44 113 272 4101

**P/N: 108978**

**Rev: C**

**08/07**

# ENGLISH

## VERSALOK™ Anchor

### Rx Only

#### Illustrations

1. End Cap
2. Threader Tab
3. Shaft
4. Pusher
5. Laser line distal to shoulder of the Pusher
6. Anchor Pin
7. Anchor Sleeve
8. Gun
9. Tension Wheel
10. Suture & Shaft Release Button

#### DESCRIPTION

The VERSALOK anchor is a two-part implant, consisting of a titanium pin and a polymer sleeve, designed to secure soft tissue to bone. When used in conjunction with #2 suture and the delivery instrumentation, the implant provides a means of firmly reattaching soft tissue to bone. The implant is supplied sterile ready to use.

The gun is a reusable instrument used to aide in the tensioning of the suture and the deployment of the anchor into bone. See IFU 108410 for cleaning and sterilization instructions for the instrument.

#### MATERIALS

**Anchor:** Polyetheretherketone (PEEK) and Titanium

**Shaft:** Stainless steel

**End Cap:** Polycarbonate

**Threader Tab:** Polycarbonate and stainless steel

1. ORTHOCORD™ Suture is a synthetic, sterile, braided, composite suture composed of dyed (D&C Violet #2 or D&C Blue #6) absorbable polydioxanone (PDS®) and un-dyed nonabsorbable polyethylene. The partially absorbable suture is coated with a copolymer of 90% caprolactone and 10% glycolide. The PDS copolymer has been shown to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

#### INDICATIONS

The VERSALOK anchor is indicated for use in the following procedures for reattachment of soft tissue to bone:

- **Shoulder:** Repair of rotator cuff tears and Biceps Tenodesis.
- **Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis, Joint Capsule Closure
- **Elbow:** Biceps Tendon Reattachment

## **CONTRAINDICATIONS**

1. Procedures other than those listed in the INDICATIONS section.
2. Pathologic conditions of bone such as cystic changes or severe osteopenia that would impair its ability to securely fix the VERSALOK anchor.
3. Pathologic changes in the soft tissues being fixated to bone that would prevent their secure fixation by the VERSALOK anchor.
4. Comminuted bone surface that would militate against secure fixation of the VERSALOK anchor.
5. Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infection, etc.
6. Conditions which tend to limit the patient's ability to restrict activities or follow directions during the healing period.
7. The VERSALOK anchor is not designed for and should never be used to attach artificial ligaments.

## **PRECAUTIONS**

1. The VERSALOK anchor is supplied STERILE, and is intended for single use only. Do not resterilize. Do not use if sterile packaging appears to be damaged.
2. The product should be stored under cool, dry conditions.
3. Inspect all instruments for damage before use. Do not attempt to repair.
4. A surgeon should not begin clinical use of the VERSALOK anchor without reviewing the Instructions for Use and practicing the procedure in a skills laboratory.
5. Difficulty tensioning sutures may occur as a result of sutures twisting around anchor, excessive suture slack around the anchor upon insertion, malleting the anchor into the bone past the laser line to the shoulder, and/or immobile Rotator Cuff.

## **WARNINGS**

1. DePuy Mitek anchors are designed to lock into cortical or cancellous bone. Bone stock must be adequate to allow proper and secure anchor

placement. Incomplete insertion or poor bone quality may result in anchor pullout.

2. Immediate range of motion should be avoided to allow biological bony/soft tissue healing.
3. This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the thoracic or lumbar spine.
4. Users should be familiar with surgical procedures and techniques involving absorbable and nonabsorbable suture before employing ORTHOCORD suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

## **INSTRUCTIONS FOR USE**

1. Remove the suture card and VERSALOK anchor from the sterile package and remove the suture from the card.
2. Pass the suture through soft tissue with the surgeon's preferred suture-passing technique.
3. Check to be sure the suture is not twisted. Load the suture onto the anchor using the attached threader tab. The threader tab pulls the suture strands through the anchor.
4. While keeping nominal tension on the suture strands, place the shaft against bone at the selected fixation site and confirm that the suture is not twisted. If the suture is twisted, untwist the suture.
5. Tap the back end of the shaft with a mallet, driving the anchor into bone. Continue tapping until contact with the distal edge of the laser line, which is directly before the shoulder of the pusher is in contact with the bone surface. In hard bone, the awl may be used first to create a starter hole. Optional: place the end cap on the shaft to increase the surface area for the mallet to strike.
6. Remove the end cap, if used, and advance the Gun (Catalog No. 214710) down the shaft until it clicks firmly into place.
7. Approximate the repair by pulling on the suture while placing gently forward pressure on the gun. The suture may be placed on the tension wheel on the VERSALOK Gun to hold the tension once it is manually set or to aide in suture tensioning. To load the suture on the tension wheel gently and equally tension the suture strands. Wrap them around the tension wheel once in a clockwise direction and then give a quick tug on the suture to set them into the tension wheel. Rotate the tension wheel clockwise to tension the suture. A maximum of two strands can be tensioned on the wheel at the same time.

8. When the desired tension is achieved confirm that the implant is still inserted or attached to the shoulder of the pusher. If the implant has backed out during the tensioning process, tap it back to the desired depth. The tension on the suture may need to be released to allow the anchor to go back in the hole.
9. Deploy the anchor by gently applying downward pressure on the handle of the VERSALOK gun and fully pulling the trigger on the VERSALOK gun.
10. If the tension wheel was used, release the suture by pressing the suture release button.
11. Slowly remove the gun from the shaft by pressing the metal shaft release button. It might be necessary to grasp the suture below the wheel while removing the gun from the shaft.
12. Place gentle tension on the suture and unthread the shaft by rotating the shaft with backwards tension in a counter-clockwise direction. Once unthreaded, completely remove the shaft and trim the excess suture.
13. Additional VERSALOK anchors may be implanted following the same technique.

## **INSTRUCTIONS FOR REINSERTING THREADER TAB**

Place the Shaft (3) and Threader Tab (2) in the same orientation as seen on the illustration. Insert each wire loop of the Threader Tab (2) through each oval hole on the anchor pin (6).

## **STORAGE**

Store in a cool dry place below 26°C (80°F), away from moisture and direct heat. Discard if open but unused.

Do not use after expiration date.

## **STERILIZATION**

The VERSALOK suture anchor is provided sterile.

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